

# Wessels company



## G-SERIES GLYMATIC SINGLE SYSTEM PACKAGE

Models: G-6 to G-18

**Submittal Sheet No.** I-1101B **Rev.** 2 2/1/2022

Job NameLocation	Submitted By Approved By	Date Date
Engineer Contractor Sales Rep.	Order No Notes	Date

#### **Description:**

The GLYMATIC is fully automated and autonomous. The unit comes complete with 110V pressurization pump with integral pressure control, pressure reducing valve and gauge, and translucent polyethylene solution container with lid designed to accommodate relief piping.

#### **Performance Limitations:**

Maximum Temperature: 130°F

Factory Discharge Pressure Setting: 12 PSIG Maximum Discharge Pressure: 60 PSIG

#### Available options:

- Alarm kit with panel and float
- Wall mounting shelf for G-6

MODEL NUMBER	PART NUMBER	TANK VOLUME (GAL)	TAGGING INFORMATION	QUANTITY
G-6	07001006	6		
G-18	07001018	18		

#### Typical Specification:

The contractor shall supply and install, as indicated on the plans and in the specifications, a prefabricated, automatic and autonomous make-up package for each glycol system. The package shall be designed to occupy a minimum amount of floor space to operate on a 110V motor, 60 Hz electrical circuit, and to maintain a fill pressure in the glycol system of \_\_\_\_\_psig. The pressurization assembly shall consist of a pressurization pump with pressure controls, a PRV, and a pressure gauge to continuously monitor PRV outlet pressure. The pressurization assembly shall be mounted on the translucent polyethylene solution container, complete with cover.

Each glycol make-up package shall be model G - \_\_\_\_\_ as manufactured by Wessels Company or approved equal.






SUBMITTAL

## **G-SERIES**GLYMATIC SINGLE SYSTEM PACKAGE

Models: G-6 to G-18

**Submittal Sheet No.** I-1101B **Rev.** 2 2/1/2022



G-6 thru G-18

### **Dimensions & Weights:**

Model Number	Dimensions in Inches				Approx
	Α	В	С	System Connection	Approx. Shipping Weight (lbs)
G-6	17 1/2	12	10	1/2 NPT	16
G-18	40		12		28